


SM Prosthetic Planning Kit

Cemented Abutment Planner, Angled Abutment Planner

SM Prosthetic Planning Kit



DIO IMPLANT

12, Centum 6 gill(1464,U-dong),
Haeundae-gu Busan, 612-020, Korea
Tel.+82-51-745-7777 Fax.+82-51-745-7778



www.dioimplant.com

DIO IMPLANT



SM Prosthetic Planning Kit

Prosthetic Planning Kit is designed to verify precision of the inserted fixture path, height of gingival, and the space between Occusal & Adjacent teeth to simulate before placing the final prosthesis.

Cemented Abutment Planner

For SM Cemented / Solid Abutment

Variety of sizes and lengths are prepared in the kit to support users to choose proper Abutments for patients.



Angled Abutment Planner

For SM Angled Abutment

Two types of angle, cuff and torx can be combined up to 24 pairs for the perfect fit.

